

Product Specification



DESCRIPTION	Potassium Bromide 98.0-100.5% Ph Eur V
-------------	--

PRODUCT CODES	GPC6945	
---------------	---------	--

EXPIRY DATE	5 Years
-------------	---------

Atom Scientific Ltd
Unit 2B
East Tame Business Park
Hyde SK14 4G
United Kingdom

T: + 44 (0) 161 3665123
E: technical@atomscientific.com

RESULTS	<p>Conforms to Ph Eur V</p> <p>Minimum Assay: >98.0-100.5% Molecular Formula: KBr Molecular Weight: 119.00 g/mol Melting Point: 730 °C Boiling Point: 1435 °C Refractive Index: 1.559 Chloride (Cl): <0.6% Sulphate (SO₄): <100ppm Iron (Fe): <20ppm Heavy Metals (as Pb): <10ppm Magnesium (Mg) and other earth metals (as Ca): <200ppm Loss on drying (3 hours): <1%</p>
---------	---

This document certifies that the above product is supplied against our approved Product Specification and complies with the company's standard specification.

A handwritten signature in black ink, appearing to read "Peter Keenan".

Peter Keenan
Commercial Director

